




## HEAD SPACE SAMPLER – Model HSS - 20

**This equipment is an accessory to be attached to and operated via Netel's Gas chromatograph instrument. The equipment contains a sample vial heater block & transfer line which can be heated up to desired temperature, automatically operated needle, a pneumatic system.**



**The equipment is connected to GC by means of carrier gas cum transfer line just ahead of column, and by Purge-In-Port, as shown in the pneumatic diagram, figure 2.**



**The carrier connections may be made either before or after the manual injector port of the GC, as indicated. The Purge connection should be made from back side for the carrier gas input. The carrier gas pressure can be adjusted with the knob behind the instrument and the pressure can be seen on front side pressure gauge. Also the carrier gas pressure input to GC can be adjusted similarly. For Netel's Gas Chromatograph the auto start electrical contact output is available for automatically starting the GC when the gas injection is being done.**



## **Technical specifications**

<b>Model Name</b>	<b>: Head Space Sampler HSS-20</b>
<b>Type</b>	<b>: Static</b>
<b>Vial Size</b>	<b>: 20 ml</b>
<b>No. of vials</b>	<b>: 20</b>
<b>Sample Loop</b>	<b>: 1 ml</b>
<b>Vial temperature range</b>	<b>: Ambient to 250°C, settable</b>
<b>Programme</b>	<b>: 1 to 20 vials, or settable</b>
<b>Valving</b>	<b>: Motorized 6 port valve</b>
<b>Valve /Loop temperature</b>	<b>: Automatically gets adjusted with Vial Temperature</b>
<b>Sample path</b>	<b>: 1/16" silylated s.s</b>
<b>Transfer Line Temperature</b>	<b>: Ambient to 250°C, Settable</b>
<b>Electronic control</b>	<b>: Parameters entry via GC front Panel</b>
<b>Physical dimensions</b>	<b>: 35.0 cm W x 31.0 cm L x 32 cm H</b>
<b>Display Characters</b>	<b>: LCD with backlit, 4 Lines x 20</b>
<b>Key Board</b>	<b>: Membrane</b>
<b>Weight</b>	<b>: 10 Kg Approx.</b>
<b>Power requirements</b>	<b>: 230 V AC, 50 Hz supply.</b>



NETEL (INDIA) LTD